

Package: sRm (via r-universe)

May 29, 2026

Type Package

Title Toolkit for SRM-MS data (MS) data

Version 0.2.6

Description A complete toolkit for working with Selective Reaction Monitoring (SRM) Mass Spectrometry (MS) data and other variants of targeted Liquid Chromatography (LC) MS data.

License GPL (>= 3)

Depends ggplot2, magrittr

Imports methods, MSnbase, tibble, purrr, dplyr, tidyr, xml2, xcms, stringr, MsCoreUtils, cli, crayon, furr, future

URL <https://github.com/wilsontom/sRm>

BugReports <https://github.com/wilsontom/sRm/issues>

LazyData TRUE

Encoding UTF-8

RoxygenNote 7.2.3

Suggests testthat, knitr, rmarkdown, patchwork

Roxygen list(markdown = TRUE)

VignetteBuilder knitr

Collate 'allGenerics.R' 'allClasses.R' 'accessor-method.R'
'adjustBaseline-method.R' 'createGroup-method.R'
'detectInstrumentModel.R' 'detectPeaks-method.R'
'fileMetaData.R' 'filterPeaks-method.R' 'freeScales.R'
'groupPeaks-method.R' 'groupSummary-method.R'
'indexRT-method.R' 'internals.R' 'keepTransitions-method.R'
'openSRM.R' 'peakAsymmetry-method.R' 'plotGroup-method.R'
'plotParent-method.R' 'plotPeakArea-method.R'
'plotSRM-method.R' 'plotSample-method.R' 'reducePeaks-method.R'
'removeSample-method.R' 'removeTransitions-method.R'
'show-method.R' 'smoothChrom-method.R'

Config/pak/sysreqs

cmake libgmp-dev make libicu-dev libuv1-dev libxml2-dev libnetcdf-dev libssl-dev zlib1g-dev

Repository <https://wilsontom.r-universe.dev>

Date/Publication 2026-05-07 20:45:36 UTC

RemoteUrl <https://github.com/wilsontom/sRm>

RemoteRef HEAD

RemoteSha d955518c522a338334588dd1d721724281a53f1c

Contents

| | |
|---------------------------------|----|
| adjustBaseline | 3 |
| createGroup | 3 |
| detectInstrumentModel | 4 |
| detectPeaks | 4 |
| fileMetaData | 5 |
| filterPeaks | 5 |
| freeScales | 6 |
| groupPeaks | 6 |
| groupSummary | 7 |
| header | 7 |
| indexRT | 8 |
| keepTransitions | 8 |
| meta | 9 |
| openSRM | 9 |
| peakAsymmetry | 10 |
| peaks | 10 |
| plotGroup | 11 |
| plotParent | 11 |
| plotPeakArea | 12 |
| plotSample | 12 |
| plotSRM | 13 |
| reducePeaks | 13 |
| removeSample | 14 |
| removeTransitions | 14 |
| show,SRM-method | 15 |
| smoothChrom | 15 |
| SRM-class | 16 |
| transitions | 16 |

Index

17

| | |
|----------------|------------------------|
| adjustBaseline | <i>Adjust Baseline</i> |
|----------------|------------------------|

Description

Adjust Baseline

Usage

```
adjustBaseline(object, hwm)
```

```
## S4 method for signature 'SRM'  
adjustBaseline(object, hwm = 6)
```

Arguments

| | |
|--------|---|
| object | a SRM object |
| hwm | a numeric value for window half width for local medians |

Value

a SRM object

| | |
|-------------|---------------------|
| createGroup | <i>Create Group</i> |
|-------------|---------------------|

Description

Create peak group based on user specified retention time

Usage

```
createGroup(object, index, rt, width, id)
```

```
## S4 method for signature 'SRM'  
createGroup(object, index, rt, width, id)
```

Arguments

| | |
|--------|--|
| object | a SRM object |
| index | a numeric value of the target transition index |
| rt | a numeric value of target Rt time (in minutes) |
| width | a numeric value of tolerated Rt width (in seconds) |
| id | a character string for peak group identification |

Value

a SRM object

detectInstrumentModel *Detect Instrument Model*

Description

Detect Instrument Model

Usage

```
detectInstrumentModel(mzml_file)
```

Arguments

mzml_file the absolute file path of a valid .mzML file

Value

a character string of the MS cvParam description for the instrument model

detectPeaks *Detect Peaks*

Description

Detect Peaks

Usage

```
detectPeaks(object, snthresh, peakwidth)
```

```
## S4 method for signature 'SRM'
```

```
detectPeaks(object, snthresh, peakwidth)
```

Arguments

object a SRM object

snthresh a numeric value for the signal-to-noise threshold to use

peakwidth a numeric vector indicating the minimum and maximum tolerated peak width

Value

a SRM object

| | |
|--------------|---|
| fileMetaData | <i>Extract meta information from .mzML file</i> |
|--------------|---|

Description

Extract a minimal amount meta data from a .mzML file

Usage

```
fileMetaData(x)
```

Arguments

x a valid .mzML file

Value

a tibble containing;

- mzML Schema
- Acquisition Date
- Acquisition Time
- Instrument Model
- File ID

| | |
|-------------|---------------------|
| filterPeaks | <i>Filter Peaks</i> |
|-------------|---------------------|

Description

Filter detected peaks based on a series of heuristic measures

Usage

```
filterPeaks(object)
```

```
## S4 method for signature 'SRM'  
filterPeaks(object)
```

Arguments

object a SRM object

Value

a SRM object

`freeScales`*Free Scales*

Description

Allow free scales in facet plots

Usage

```
freeScales(ggplot)
```

Arguments

`ggplot` a ggplot object

Value

a ggplot object

`groupPeaks`*Group Peaks*

Description

Group detected peaks based on retention time. Based on `MsCoreUtils::group`

Usage

```
groupPeaks(object, rt_tolerance)
```

```
## S4 method for signature 'SRM'  
groupPeaks(object, rt_tolerance)
```

Arguments

`object` a SRM object

`rt_tolerance` a numeric value for the tolerate retention grouping

Value

a SRM object

| | |
|---------------------------|----------------------|
| <code>groupSummary</code> | <i>Group Summary</i> |
|---------------------------|----------------------|

Description

Create a summary table for peak groups

Usage

```
groupSummary(object)
```

```
## S4 method for signature 'SRM'  
groupSummary(object)
```

Arguments

`object` a SRM object

Value

a tibble of group summary

| | |
|---------------------|---------------|
| <code>header</code> | <i>header</i> |
|---------------------|---------------|

Description

List object header information

Usage

```
header(object)
```

```
## S4 method for signature 'SRM'  
header(object)
```

Arguments

`object` a SRM S4 object

| | |
|---------|---|
| indexRT | <i>Find Target Retention Time (Rt) for a transition</i> |
|---------|---|

Description

Find Target Retention Time (Rt) for a transition

Usage

```
indexRT(object, index, sampleName, n)
```

```
## S4 method for signature 'SRM'
```

```
indexRT(object, index, sampleName, n = 1)
```

Arguments

| | |
|------------|---|
| object | a SRM object |
| index | a numeric value for the transition index |
| sampleName | a character string |
| n | a numeric value for the number of Rts to return (if > 1 is available) |

Value

a numeric vector or Rt targets

| | |
|-----------------|-------------------------|
| keepTransitions | <i>Keep Transitions</i> |
|-----------------|-------------------------|

Description

Keep Transitions

Usage

```
keepTransitions(object, index_keep)
```

```
## S4 method for signature 'SRM'
```

```
keepTransitions(object, index_keep = c())
```

Arguments

| | |
|------------|--|
| object | a SRM object |
| index_keep | a numeric vector of transition index to retain |

Value

a SRM object

meta

meta

Description

List sample meta data

Usage

```
meta(object)
```

```
## S4 method for signature 'SRM'  
meta(object)
```

Arguments

object a SRM S4 object

openSRM

Open SRM Files

Description

Open and parse SRM files into an S4 SRM Object

Usage

```
openSRM(files)
```

Arguments

files a character vector of absolute file paths of SRM files in .mzML format

Value

an SRM object

| | |
|---------------|-----------------------|
| peakAsymmetry | <i>Peak Asymmetry</i> |
|---------------|-----------------------|

Description

Peak Asymmetry

Usage

```
peakAsymmetry(object)
```

```
## S4 method for signature 'SRM'  
peakAsymmetry(object)
```

Arguments

object a SRM object

Value

a SRM object

| | |
|-------|--------------|
| peaks | <i>peaks</i> |
|-------|--------------|

Description

List detected peaks

Usage

```
peaks(object)
```

```
## S4 method for signature 'SRM'  
peaks(object)
```

Arguments

object a SRM S4 object

| | |
|-----------|-------------------|
| plotGroup | <i>Plot Group</i> |
|-----------|-------------------|

Description

Plot Group

Usage

```
plotGroup(object, group)

## S4 method for signature 'SRM'
plotGroup(object, group)
```

Arguments

| | |
|--------|-----------------------------------|
| object | a SRM object |
| group | a character string of the GroupID |

Value

a ggplot plot object

| | |
|------------|--------------------|
| plotParent | <i>Plot Parent</i> |
|------------|--------------------|

Description

Plot Parent

Usage

```
plotParent(object, parentMass)

## S4 method for signature 'SRM'
plotParent(object, parentMass)
```

Arguments

| | |
|------------|--|
| object | a SRM object |
| parentMass | a numeric value of the parent mass (Q1) to extract |

Value

a ggplot plot object

plotPeakArea *Plot Peak Area*

Description

Plot Peak Area

Usage

```
plotPeakArea(object, index, sampleName)
```

```
## S4 method for signature 'SRM'
```

```
plotPeakArea(object, index, sampleName)
```

Arguments

| | |
|------------|---|
| object | a SRM object |
| index | a numeric value of the transition index to plot |
| sampleName | a character string of sampleName to plot |

Value

a ggplot plot object

plotSample *Plot Sample*

Description

Plot Sample

Usage

```
plotSample(object, sampleName, polarity)
```

```
## S4 method for signature 'SRM'
```

```
plotSample(object, sampleName, polarity)
```

Arguments

| | |
|------------|--|
| object | a SRM object |
| sampleName | a character string of a valid sampleName |
| polarity | a character string of either pos or neg for positive and negative ionisation mode respectively |

Value

a ggplot plot object

| | |
|---------|-----------------|
| plotSRM | <i>Plot SRM</i> |
|---------|-----------------|

Description

Plot SRM

Usage

```
plotSRM(object, index, type = "overlay")
```

```
## S4 method for signature 'SRM'
plotSRM(object, index, type = "overlay")
```

Arguments

| | |
|--------|---|
| object | a SRM object |
| index | a numeric value of the transition index to plot |
| type | a character string of either overlay or facet |

Value

a ggplot plot object

| | |
|-------------|---------------------|
| reducePeaks | <i>Reduce Peaks</i> |
|-------------|---------------------|

Description

Reduce detected peaks using a relative area threshold

Usage

```
reducePeaks(object, relative_area)
```

```
## S4 method for signature 'SRM'
reducePeaks(object, relative_area = 0.05)
```

Arguments

| | |
|---------------|--|
| object | a SRM object |
| relative_area | a numeric value between 0 and 1, for the threshold of relative peak area |

Value

a SRM object

| | |
|--------------|----------------------|
| removeSample | <i>Remove Sample</i> |
|--------------|----------------------|

Description

Remove Sample

Usage

```
removeSample(object, sampleName)
```

```
## S4 method for signature 'SRM'
removeSample(object, sampleName)
```

Arguments

| | |
|------------|--|
| object | a SRM object |
| sampleName | a character string of a sampleName to remove |

Value

a SRM object

| | |
|-------------------|-------------------------|
| removeTransitions | <i>Keep Transitions</i> |
|-------------------|-------------------------|

Description

Keep Transitions

Usage

```
removeTransitions(object, index_out)
```

```
## S4 method for signature 'SRM'
removeTransitions(object, index_out = c())
```

Arguments

| | |
|------------|--|
| object | a SRM object |
| index_keep | a numeric vector of transition index to remove |

Value

a SRM object

| | |
|------------------|-----------------|
| show, SRM-method | <i>show-SRM</i> |
|------------------|-----------------|

Description

show-SRM

Usage

```
## S4 method for signature 'SRM'
show(object)
```

Arguments

object a SRM object

| | |
|-------------|-----------------------------|
| smoothChrom | <i>Smooth Chromatograms</i> |
|-------------|-----------------------------|

Description

Smooth Chromatograms

Usage

```
smoothChrom(object, cf)

## S4 method for signature 'SRM'
smoothChrom(object, cf)
```

Arguments

object a SRM object

cf matrix, a coefficient matrix generated by MsCoreUtils::coefMA or MsCore-
Utils::coefWMA

Value

a SRM object

 SRM-class

SRM

Description

A S4 class to store single reaction monitoring (SRM) mass spectrometry (MS) data

Slots

meta a tibble of a minimal amount of meta data.

transitions a tibble of all transitions names

chroms a tibble of retention time (rt) and intensity (int) values for all imported SRM chromatograms

peaks a tibble of detected peaks

groups a tibble of peak information for user defined groups

header a tibble of formatted header data from mzR: :header()

 transitions

Accessor Methods for the SRM S4 class

Description

List all transitions

Usage

```
transitions(object)
```

```
## S4 method for signature 'SRM'
transitions(object)
```

Arguments

object a SRM S4 object

Index

[adjustBaseline](#), [3](#)
[adjustBaseline](#), SRM-method
 ([adjustBaseline](#)), [3](#)

[createGroup](#), [3](#)
[createGroup](#), SRM-method ([createGroup](#)), [3](#)

[detectInstrumentModel](#), [4](#)
[detectPeaks](#), [4](#)
[detectPeaks](#), SRM-method ([detectPeaks](#)), [4](#)

[fileMetaData](#), [5](#)
[filterPeaks](#), [5](#)
[filterPeaks](#), SRM-method ([filterPeaks](#)), [5](#)
[freeScales](#), [6](#)

[groupPeaks](#), [6](#)
[groupPeaks](#), SRM-method ([groupPeaks](#)), [6](#)
[groupSummary](#), [7](#)
[groupSummary](#), SRM-method ([groupSummary](#)),
 [7](#)

[header](#), [7](#)
[header](#), SRM-method ([header](#)), [7](#)

[indexRT](#), [8](#)
[indexRT](#), SRM-method ([indexRT](#)), [8](#)

[keepTransitions](#), [8](#)
[keepTransitions](#), SRM-method
 ([keepTransitions](#)), [8](#)

[meta](#), [9](#)
[meta](#), SRM-method ([meta](#)), [9](#)

[openSRM](#), [9](#)

[peakAsymmetry](#), [10](#)
[peakAsymmetry](#), SRM-method
 ([peakAsymmetry](#)), [10](#)
[peaks](#), [10](#)

[peaks](#), SRM-method ([peaks](#)), [10](#)
[plotGroup](#), [11](#)
[plotGroup](#), SRM-method ([plotGroup](#)), [11](#)
[plotParent](#), [11](#)
[plotParent](#), SRM-method ([plotParent](#)), [11](#)
[plotPeakArea](#), [12](#)
[plotPeakArea](#), SRM-method ([plotPeakArea](#)),
 [12](#)
[plotSample](#), [12](#)
[plotSample](#), SRM-method ([plotSample](#)), [12](#)
[plotSRM](#), [13](#)
[plotSRM](#), SRM-method ([plotSRM](#)), [13](#)

[reducePeaks](#), [13](#)
[reducePeaks](#), SRM-method ([reducePeaks](#)), [13](#)
[removeSample](#), [14](#)
[removeSample](#), SRM-method ([removeSample](#)),
 [14](#)
[removeTransitions](#), [14](#)
[removeTransitions](#), SRM-method
 ([removeTransitions](#)), [14](#)

[show](#), SRM-method, [15](#)
[smoothChrom](#), [15](#)
[smoothChrom](#), SRM-method ([smoothChrom](#)), [15](#)
[SRM-class](#), [16](#)

[transitions](#), [16](#)
[transitions](#), SRM-method ([transitions](#)), [16](#)